

# PathoCatch<sup>®</sup>

## HBsAg Test Kit

Behind every diagnosis, There is a proven performance.



### Simplified, Standardized & Rapid testing solutions.

Hepatitis B caused by Hepatitis B Virus (HBV) is responsible for chronic infection of the liver and increases the chances of cirrhosis and liver cancer. As per the estimates, there were nearly 296 million HBV cases, and resulting in the death toll of about 0.82 million in 2019. Around 1.5 million new infections are estimated to originate each year. Early detection & timely testing are critical for the management of Hepatitis B infection.

Source: WHO

## PathoCatch® HBsAg Test Kit

Early screening has positively affected the risk and transmission of HBV infection. In its initial stage, early detection assists in managing the infection. As per the recommendations, rapid antigen testing for HBsAg is the method of choice for accurate detection of HBV. The test functions by identifying surface antigens in the bloodstream.

PathoCatch® HBsAg test is a rapid antigen test kit that accurately detects Hepatitis B Surface antigen in whole blood, serum and plasma. It is a simple, easy-to-use test kit. Mylab's HBV detection kits are accurate, sensitive, and highly specific for identifying HBV thereby assisting in eliminating the threat of HBV infection.

### Kit Features



#### Accurate results

Highly accurate and reliable results



#### Simple sampling

Sample: whole blood, serum or plasma

Minimal sample of 50 µl



#### Easy to use

Simple single-step process

Result in 20 mins

Easy to interpret results



#### Stable to store

Shelf life of 24 months.

Stable at ambient temperature

### Ordering Information

Product Name	Cat. No.	Reaction Size/Kit
PathoCatch® HBsAg Test	2013A025	25 Tests
PathoCatch® HBsAg Test	2013A050	50 Tests
PathoCatch® HBsAg Test	2013A100	100 Tests



### Related Products for HBs Ag Test Kit



**Mylab Mybox+™**  
Dry Luminescence Assay Analyzer



**Pathodetect®**  
Automated HBV PCR Kit



**Compact DX®**  
Automated sample to report system



#### Mylab Discovery Solutions Pvt. Ltd.

Global Innovation Center (GIC) 5th Floor, Amar Paradigm, Survey No - 110/11/3, Baner Rd, Pune - 411045, Maharashtra, India. | Email: info@mylabglobal.com | Web: www.mylabglobal.com

PathoDetect | PathoCatch | NATSpert | CompactDX | CompactXL | OncoScreen | DiscoverSeries | AgroDex | VetScreen | CoviSelf | CoviSwift

**Disclaimer:** All rights reserved. All content (texts, trademarks, illustrations, photos, graphics, files, designs, arrangements, instruments etc.) are protected by our Intellectual Property (IP) including Patents, Copyright, Trademark and other protective laws. Any attempt for infringement of IP vested in our products, for any purpose whatsoever, without taking prior written license from Mylab would amount to infringement and may be harmful involving high risk to life and may result in fatal injury. For enquiries related to IP licensing, please write to legal@mylabglobal.com