

PathoCatch[®]

HCV Ab Test Kit

Behind every diagnosis, There is a proven performance.



Simplified, Standardized & Rapid testing solutions.

Hepatitis C Virus (HCV) is a type of hepatotropic virus causing acute and chronic liver infection. As per the WHO data, 58 million individuals were suffering from chronic HCV infection in 2019, with more than 0.29 million deaths yearly. Although the global eradication of HCV has been in overdrive, there are still gaps in the testing of HCV.

Source: WHO

PathoCatch® HCV Ab Test Kit

Early detection enables immediate management of the infection in its initial stage. Screening has massive implications for reducing the risk of HCV infection and transmission to other people. It starts with antibody testing using rapid antibody tests. They function by detecting antibodies to HCV in the bloodstream.

Mylab with its rapid antibody test kits provides unprecedented and accurate testing results for the detection of HCV infection. Mylab with its specialized diagnostic services has been a major player in the diagnostic space. The sincere effort in reducing the risk of deadly infections and providing affordable diagnostic solutions has been a major outlook of Mylab.

Kit Features



Accurate results

Detection of all the subtypes of HCV

Highly accurate and reliable results



Simple sampling

Sample: whole blood, serum or plasma samples

Minimal sample of 25 µl



Easy to use

Simple two-step process

Result in 20 mins

Easy to interpret results



Stable to store

Shelf life of 24 months.

Stable at ambient temperature

Ordering Information

Product Name	Cat. No.	Reaction Size/Kit
PathoCatch® HCV Ab Test	2014A025	25 Tests
PathoCatch® HCV Ab Test	2014A050	50 Tests
PathoCatch® HCV Ab Test	2014A0100	100 Tests



Related Products for HCV Ab Test Kit



Mylab Mybox+™
Dry Luminescence Assay Analyzer



Pathodetect®
automated HCV PCR Kit



Compact DX®
Automated sample to report system



Mylab Discovery Solutions Pvt. Ltd.

Global Innovation Center (GIC) 5th Floor, Amar Paradigm, Survey No - 110/11/3, Baner Rd, Pune - 411045, Maharashtra, India.

Email: info@mylabglobal.com | **Web:** www.mylabglobal.com

PathoDetect | PathoCatch | NATSpert | CompactDX | CompactXL | OncoScreen | DiscoverSeries | AgroDex | VetScreen | CoviSelf | CoviSwift

Disclaimer: All rights reserved. All content (texts, trademarks, illustrations, photos, graphics, files, designs, arrangements, instruments etc.) are protected by our Intellectual Property (IP) including Patents, Copyright, Trademark and other protective laws. Any attempt for infringement of IP vested in our products, for any purpose whatsoever, without taking prior written license from Mylab would amount to infringement and may be harmful involving high risk to life and may result in fatal injury. For enquiries related to IP licensing, please write to legal@mylabglobal.com