

PathoCatch[®]

HIV 1/2 Ab Test Kit

Behind every diagnosis, There is a proven performance.



Simplified, Standardized & Rapid testing solutions.

Human Immunodeficiency Virus (HIV) belongs to the group of lentiviruses capable of disrupting the normal functioning of the immune system, The advanced stage of HIV is usually characterized by acquired immunodeficiency syndrome (AIDS). As of 2022, more than 39 million people are living with HIV, and around 0.63 million individuals have died. Hence, early detection and diagnosis are one of the few methods for the prevention of HIV infection.

Source: WHO

PathoCatch® HIV 1/2 Ab Test Kit

HIV testing is one of the crucial aspects of the detection and management of HIV infection. It consists of laboratory diagnosis via Rapid Diagnostic Tests (RDTs). The test detects the presence of HIV-1 and HIV-2 antibodies. RDTs are highly accurate, reliable, efficient, and easy to perform.

Mylab provides state-of-the-art rapid diagnostic kits for the detection of HIV. The company delves into specialized, reliable and high-quality diagnostic solutions. The impeccable accuracy and specificity for the detection of pathogens are the unique features of Mylab's detection kits.

Kit Features



Accurate results

Detection of HIV 1 and HIV 2

Highly accurate & reliable results



Simple sampling

Sample: whole blood, serum or plasma samples

Minimal sample of 25 µl



Easy to use

Simple two-step process

Result in 20 mins

Easy to interpret results



Stable to store

Shelf life of 24 months.

stable at ambient temperature

Ordering Information

Product Name	Cat. No.	Reaction Size/Kit
PathoCatch® HIV1/2 Ab Test	2015A025	25 Tests
PathoCatch® HIV1/2 Ab Test	2015A050	50 Tests
PathoCatch® HIV1/2 Ab Test	2015A0100	100 Tests



Related Products for HIV 1/2 Ab Test Kit



Mylab Mybox+™
Luminescence Assay Analyzer



Pathodetect®
Automated HIV PCR Kit



Compact DX®
Automated sample to report system



Mylab Discovery Solutions Pvt. Ltd.

Global Innovation Center (GIC), 5th Floor, Amar Paradigm, Survey No - 110/11/3, Baner Road, Pune - 411045, Maharashtra, INDIA. | Email: info@mylabglobal.com | Web: www.mylabglobal.com

PathoDetect | PathoCatch | NATSpert | CompactDX | CompactXL | OncoScreen | DiscoverSeries | AgroDex | VetScreen | CoviSelf | CoviSwift

Disclaimer: All rights reserved. All content (texts, trademarks, illustrations, photos, graphics, files, designs, arrangements, instruments etc.) are protected by our Intellectual Property (IP) including Patents, Copyright, Trademark and other protective laws. Any attempt for infringement of IP vested in our products, for any purpose whatsoever, without taking prior written license from Mylab would amount to infringement and may be harmful involving high risk to life and may result in fatal injury. For enquiries related to IP licensing, please write to legal@mylabglobal.com