

# PathoCatch™ COVID-19 Ag Detection Kit



COVID-19 is an acute respiratory infectious disease. Currently, the patients are infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. Mylab has developed COVID-19 Ag Card test which is used for in vitro qualitative detection of the antigen of novel Coronavirus in human throat swabs or nasal swabs. It employs immunochromatography for detection.

## Product Features:



**Easier:** No special equipment needed, Intuitive visual interpretation.



**Rapid:** Quick sampling, Result in 20 minutes.



**Accurate:** Results validated using PCR

**Application:** For suspicious patients with symptoms, mild symptoms, or even without symptoms, also for testing people with close contact of infected patients and people under quarantine control.

**Specimen:** Nasopharyngeal swab

All necessary reagents provided & no additional equipment needed

### Mylab Discovery Solutions Pvt. Ltd.

Plot No. 99-B, Lonavala, Industrial Co-operative Estate Ltd.,  
Nangargaon, Lonavala, Pune, Maharashtra 410401, INDIA.

**Email:** info@mylabglobal.com | **Web:** www.mylabglobal.com

### ORDER INFORMATION

Pack Size	REF
25 Tests	CVGC 25
50Tests	CVGC 25